

Automate acid/base titrations for simple accuracy

Potentiometric titrator for accurate pH titrations with easily reportable results.



Markets

- Food and Beverage
- Environmental
- Water
- Wastewater
- Industrial
- Consumer Products
- Academic
- Research

The Thermo Scientific™ Orion™ T910 pH titrators are designed to increase your laboratory productivity by automating titrations. Our core electrochemistry technology is integrated with a state-of-the-art reagent dispensing system to create modern, simplified automated titrators designed to make performing titrations easier, more reliable and more reproducible than manual titrations.

Use the Orion Star T910 pH titrator for dedicated acid-base titrations including:

- Titratable acidity of juices and wines
- Acidity of food products
- Alkalinity of waters
- Acidity and alkalinity of consumer products
- Total acid number (TAN)
- Total base number (TBN)

Titration techniques include equivalence point titrations and preset pH endpoint titrations for versatile sample analysis. The automated titrant delivery system with high-accuracy burette and dispenser provides consistent results. Automatic results calculations simplifies analysis and reduces computational errors.

This compact titrator is designed to:

- Offer easy-to-use navigation from setup to live titration analysis to data transfer
- Showcase titrations, onscreen instructions and help menus on a large color graphic touchscreen display
- Create and save up to ten user-defined methods
- Save up to 100 titration data sets with date/time stamp
- Facilitate data transfer to printer, computer or USB drive

Durable Titrators with Rugged Design

Benefits of Orion Star T910 pH Titrator:

- **Accurate, Reproducible Results:** Automating the titration workflow enables analysis to be run using the same setup parameters and calculations every time.
- **Fast Setup Time:** Setup and save electrode, titrant and titration parameters once using methods.
- **Prioritize Lab Throughput:** Start the titration, and let the instrument run unattended.
- **Minimal Cost and Downtime:** Easily replace titrator parts to minimize maintenance downtime while long lasting electrodes may reduce the cost of ownership.
- **Securely Document Test Results:** View results live, in the summary table and within the 100-point data log.
- **Applications Database:** Tested and validated preprogrammed methods can be easily loaded onto your titrator via USB flash drive.



Titrator connections have clearly labeled functions, including three USB ports for printer, USB flash drive, and computer software connections.



Minimize downtime with easy-to-install and easy-to-replace burettes, tubing, dispenser and additional high-wear components.




Electrode holder provides secure, smooth electrode movement throughout the measurement process.



Pair with premium ROSS pH electrode technology to provide a complete measurement system optimized for accurate, reproducible sample analysis.



Bottle holder helps secure titrant to prevent spills and can be easily disconnected from the titrator for flexible placement on the lab bench.



Color-coded tubing and connectors simplify tubing installation and maintenance.



Stirrer probe provides consistent sample mixing, helps prevent damage to electrodes and is easy to rinse between samples.

Product Specifications

Orion Star T910 pH Titrator

Titration Technique	Equivalence point or preset endpoint
Equivalence Points	1 or 2
Preset Endpoints	1, 2 or 3
Titration Types	Direct or back titration
Blank Options	Fixed value or variable value using titration
Cycles Per Titration Run	Up to 5 cycles with option to exclude cycles from average results and RSD calculations
Titrant Determination	Standardization titration or manual concentration input
Titration Process Control	Routine, quick, careful or user defined
Titration Precision	±0.5 % RSD, dependent on environment and handling conditions
Stirrer Probe Speeds	5 selectable speeds, 250 to 3700 RPM
Sample ID	Auto-incremental, manual or off
Titration Setup Wizard	Yes, using Start a New Titration workflow
Methods	Up to 10 customizable methods with optional password protection
Methods Transfer	Import/export via USB flash drive, summary to computer software or compact printer
Data Logs	Sample titration, titrant standardization, calibration, direct measure – 100 data sets each
Data Log Export	CSV or Report (PDF) file, short or long format
Time and Date	Yes, with nonvolatile battery backup
Direct Measure Mode	pH with mV and temperature
pH Range / Resolution / Relative Accuracy	-2.000 to 20.000 pH / 0.001, 0.01, 0.1 (user selectable) / ±0.002 pH
pH Calibration Mode	1 to 5 point
mV Range / Resolution / Relative Accuracy	-2000.0 to +2000.0 mV / 0.1 mV / ±0.2 mV
Temperature Range / Resolution / Relative Accuracy	-5.0 to 100.0 °C; 23.0 to 212 °F / 0.1 °C; 0.1 °F / ±0.2 °C
Temperature Input	Manual or automatic with optional 1 point ATC probe offset calibration
Display Type	5.7" color capacitive touchscreen, 640 x 480 resolution, with lab glove compatibility
Display Backlight	Yes, adjustable
Languages	Chinese, English, French, German, Italian, Japanese, Korean, Portuguese, Spanish
Titration Setup Wizard	Yes
Audible Notifications	Titration Cycle Complete, Maximum Titrant Volume, Data Log Full, Calibration Due, Maintenance Reminder
Updateable Firmware	Yes, using USB flash drive
Burette Sizes	10 mL, 20 mL (included), 50 mL
Burette Resolution	Advanced micro-step technology provides 25,600 micro-steps per motor revolution for smooth and precise burette positioning (2 million micro-steps over full burette stroke range)
Burette Dosing Accuracy	Meets ISO8655-3 requirements
Burette Functions	Automatic rinse cycles for flushing and discrete dispense volume with continuous option
Certifications	CE, TUV 3-in-1, FCC, EN/EIC61326-1, IEC 61010, IP-51
Dimensions	10" x 16" x 14" (L x W x H); 25.4 cm x 40.6 cm x 35.6 cm (L x W x H)
Weight	12.5 lbs; 5.67 kg
Power Requirements	100 to 240 V, 50/60 Hz
Warranty	1 year





Expand the functionality of your Orion Star T910 pH titrator with key accessories designed to help increase your lab productivity. Easy to install replacement parts minimize maintenance downtime and ensure your Orion Star T910 pH titrator is ready for sample analysis.

- For a bundled discount, pair the Orion Star T910 pH titrator with any premium ROSS pH electrodes or select a preconfigured titrator and electrode kit.
- The 8102BNUWP ROSS Ultra pH electrode and 327007MD stainless steel ATC temperature probe are ideal for most standard pH titrations and preset endpoint titrations including temperature.
- The 8172BNWP ROSS Sure-Flow electrode is designed for challenging samples that can clog standard pH electrodes or that require a fast-flowing sleeve-type junction.

Ordering Information	
Cat. No.	Description
STARA-106	Orion compact ink ribbon printer, 100-240 V
START-PC1	Orion Star T900 series computer software
START-B10	Orion Star T900 series 10 mL burette
START-B20	Orion Star T900 series 20 mL burette
START-B50	Orion Star T900 series 50 mL burette
START-TB1	Orion Star T900 series standard tubing kit
START-TB2	Orion Star T900 series light blocking tubing kit
START-TB4	Orion Star T900 series burette tubing
START-TB5	Orion Star T900 series reagent bottle tubing
START-TB6	Orion Star T900 series dispenser tubing
START-BT1	Orion Star T900 series 1 L plastic bottle
START-BT2	Orion Star T900 series 1 L glass amber bottle
START-BT3	Orion Star T900 series 1 L plastic bottles, 12 pack
START-BT4	Orion Star T900 series 1 L glass amber bottles, 12 pack
START-CP1	Orion Star T900 series GL38 reagent bottle cap
START-CP2	Orion Star T900 series GL45 reagent bottle cap
START-CP3	Orion Star T900 series Orion pint bottle cap
START-DS1	Orion Star T900 series dispenser probe
START-DVK	Orion Star T900 series dispense verification kit
START-B00	Orion Star T900 series burette cover
START-BT0	Orion Star T900 series reagent bottle holder
START-EH1	Orion Star T900 series electrode holder assembly
START-EH2	Orion Star T900 series stopper for electrode holder
START-EH3	Orion Star T900 series cable management accessory
START-PS1	Orion Star T900 series 110 V US/Japan power supply
START-PS2	Orion Star T900 series 220 V Euro power supply
START-PS3	Orion Star T900 series 240 V UK/Singapore power supply
START-PS4	Orion Star T900 series 230 V Australia/New Zealand power supply
START-PS5	Orion Star T900 series 220 V China power supply
START-TB3	Orion Star T900 series drying tube
START-UM1	Orion Star T900 series user manual on USB
START-USB	Orion Star T900 series USB computer cable
START-VK1	Orion Star T900 series valve replacement kit

Ordering Information	
Cat. No.	Description
START9100	Orion Star T910 pH titrator, includes 20 mL burette, stirrer probe, dispenser probe, standard tubing kit, 1 L plastic bottle, GL38 bottle cap with desiccant tube, computer cable, literature on USB drive, 110-240 V power supplies
START9101	Orion Star T910 pH titrator standard ROSS kit, includes 8102BNUWP ROSS Ultra pH electrode, 927007MD ATC probe, 20 mL burette, stirrer probe, dispenser probe, standard tubing kit, 1 L plastic bottle, GL38 bottle cap with desiccant tube, computer cable, literature on USB drive, 110-240 V power supplies
START9102	Orion Star T910 pH titrator Sure-Flow ROSS kit, includes 8172BNWP ROSS Sure-Flow pH electrode, 927007MD ATC probe, 20 mL burette, stirrer probe, dispenser probe, standard tubing kit, 1 L plastic bottle, GL38 bottle cap with desiccant tube, computer cable, literature on USB drive, 110-240 V power supplies

Water and Lab Products

Australia: (613) 9757-4300 In Australia: (1300) 735-295 **China:** (86) 21-6865-4588
Germany: (49) 6184-90-6321 **India:** (91) 22-6716-2261/2247 **Japan:** (81) 045-453-9175
North America: 1-978-232-6000 Toll Free: 1-800-225-1480 **Singapore:** (65) 6778-6876



Tel: 45 95 04 10
 info@buhl-bonsoe.dk
 www.buhl-bonsoe.dk

Find out more at thermofisher.com/T900titratorseries

